

FUCA1 Polyclonal Antibody PE Conjugated

Catalog No: #C91858PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FUCA1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	321-420/466
Conjugates	PE
Target Name	FUCA1
Other Names	FUCA1
Accession No.	Uniprot:P04066Gene ID:2517
Uniprot	P04066
GeneID	2517
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only