

C3a Receptor Polyclonal Antibody Biotin Conjugated

Catalog No: #C91866Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C3a Receptor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-482/482
Conjugates	Bio
Target Name	C3a Receptor
Other Names	C3a Receptor
Accession No.	Uniprot:Q16581Gene ID:719
Uniprot	Q16581
GeneID	719
Calculated MW	54
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only