

Noggin Polyclonal Antibody Biotin Conjugated

Catalog No: #C91871Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Noggin Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 28-130/232 |
| Conjugates | Bio |
| Target Name | Noggin |
| Other Names | Noggin |
| Accession No. | Uniprot:Q13253Gene ID:9241 |
| Uniprot | Q13253 |
| GeneID | 9241 |
| Calculated MW | 23 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only