

GGT1 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C91872PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | GGT1 Polyclonal Antibody PE-Cy5 Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 351-460/569 |
| Conjugates | PE-Cy5 |
| Target Name | GGT1 |
| Other Names | GGT1 |
| Accession No. | Uniprot:P19440Gene ID:2678 |
| Uniprot | P19440 |
| GeneID | 2678 |
| Calculated MW | 61 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Note: This product is for in vitro research use only