

## PCGF4/BMI1 Polyclonal Antibody FITC Conjugated

Catalog No: #C91878FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PCGF4/BMI1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	61-160/326
Conjugates	FITC
Target Name	PCGF4/BMI1
Other Names	PCGF4/BMI1
Accession No.	Uniprot:P35226Gene ID:648
Uniprot	P35226
GeneID	648
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1ug/Test IF=1:200-800

Note: This product is for in vitro research use only