

# Phospho-HDAC4 (Ser246) / HDAC5 (Ser259) / HDAC7 (Ser155) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91887Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospho-HDAC4 (Ser246) / HDAC5 (Ser259) / HDAC7 (Ser155) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	TA(p-S)EP
Conjugates	Bio
Target Name	Phospho-HDAC4 (Ser246) / HDAC5 (Ser259) / HDAC7 (Ser155)
Other Names	Phospho-HDAC4 (Ser246) / HDAC5 (Ser259) / HDAC7 (Ser155)
Accession No.	Uniprot:P56524Gene ID:9759
Uniprot	P56524
GeneID	9759
Calculated MW	119
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only