

Phospho-Beta-Catenin (Ser33 / Ser37) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91890Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-Beta-Catenin (Ser33 / Ser37) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	LD(p-S)GIH(p-S)GA
Conjugates	Bio
Target Name	Phospho-Beta-Catenin (Ser33 / Ser37)
Other Names	Phospho-Beta-Catenin (Ser33 / Ser37)
Accession No.	Uniprot:P35222Gene ID:1499
Uniprot	P35222
GeneID	1499
Calculated MW	86
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only