

phospho-Cortactin (Thr466) Polyclonal Antibody FITC Conjugated



Catalog No: #C91892FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-Cortactin (Thr466) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu
Immunogen Type	YA(p-T)EA
Conjugates	FITC
Target Name	phospho-Cortactin (Thr466)
Other Names	phospho-Cortactin (Thr466)
Accession No.	Uniprot:Q14247Gene ID:2017
Uniprot	Q14247
GeneID	2017
Calculated MW	61
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only