

Phospho-Doublecortin (Ser128) Polyclonal Antibody PE Conjugated



Catalog No: #C91893PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-Doublecortin (Ser128) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	AL(p-S)NE
Conjugates	PE
Target Name	Phospho-Doublecortin (Ser128)
Other Names	Phospho-Doublecortin (Ser128)
Accession No.	Uniprot:O43602Gene ID:1641
Uniprot	O43602
GeneID	1641
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only