

Phospho-FOXO1 / FOXO3 / FOXO4 (phospho T24 / T32) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91898Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-FOXO1 / FOXO3 / FOXO4 (phospho T24 / T32) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	SC(p-T)WP
Conjugates	Bio
Target Name	Phospho-FOXO1 / FOXO3 / FOXO4 (phospho T24 / T32)
Other Names	Phospho-FOXO1 / FOXO3 / FOXO4 (phospho T24 / T32)
Accession No.	Uniprot:Q12778Gene ID:2308
Uniprot	Q12778
GeneID	2308
Calculated MW	74
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only