

# Phospho-FoxO1 (Ser256) Polyclonal Antibody PE Conjugated



Catalog No: #C91900PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospho-FoxO1 (Ser256) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	AA(p-S)MD
Conjugates	PE
Target Name	Phospho-FoxO1 (Ser256)
Other Names	Phospho-FoxO1 (Ser256)
Accession No.	Uniprot:Q12778Gene ID:2308
Uniprot	Q12778
GeneID	2308
Calculated MW	72
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only