

phospho-FOXO4 (Ser262) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91901Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-FOXO4 (Ser262) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	SS(p-S)NA
Conjugates	Bio
Target Name	phospho-FOXO4 (Ser262)
Other Names	phospho-FOXO4 (Ser262)
Accession No.	Uniprot:P98177Gene ID:4303
Uniprot	P98177
GeneID	4303
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only