

Phospho-GCN2 (Thr899) Polyclonal Antibody PE Conjugated

Catalog No: #C91905PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-GCN2 (Thr899) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	HL(p-T)GM
Conjugates	PE
Target Name	Phospho-GCN2 (Thr899)
Other Names	Phospho-GCN2 (Thr899)
Accession No.	Uniprot:Q9P2K8Gene ID:440275
Uniprot	Q9P2K8
GeneID	440275
Calculated MW	189
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test

Note: This product is for in vitro research use only