

# Phospho-IRAK1 (Ser376) Polyclonal Antibody PE Conjugated

Catalog No: #C91914PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospho-IRAK1 (Ser376) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	QS(p-S)TV
Conjugates	PE
Target Name	Phospho-IRAK1 (Ser376)
Other Names	Phospho-IRAK1 (Ser376)
Accession No.	Uniprot:Q62406Gene ID:16179
Uniprot	Q62406
GeneID	16179
Calculated MW	78
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only