

Phospho-IRF3 (Ser396) Polyclonal Antibody PE Conjugated

Catalog No: #C91917PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-IRF3 (Ser396) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	LHI(p-S)NS
Conjugates	PE
Target Name	Phospho-IRF3 (Ser396)
Other Names	Phospho-IRF3 (Ser396)
Accession No.	Uniprot:Q14653Gene ID:3661
Uniprot	Q14653
GeneID	3661
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only