

Phospho-IRF7 (Ser471 / Ser472) Polyclonal Antibody PE Conjugated



Catalog No: #C91918PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-IRF7 (Ser471 / Ser472) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	GV(p-S)(p-S)LD
Conjugates	PE
Target Name	Phospho-IRF7 (Ser471 / Ser472)
Other Names	Phospho-IRF7 (Ser471 / Ser472)
Accession No.	Uniprot:Q92985Gene ID:3665
Uniprot	Q92985
GeneID	3665
Calculated MW	54
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only