

# Phospho-Jak1 (Tyr1034 / Tyr1035) Polyclonal Antibody PE Conjugated

Catalog No: #C91920PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospho-Jak1 (Tyr1034 / Tyr1035) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	KE(p-Y)(p-Y)TV
Conjugates	PE
Target Name	Phospho-Jak1 (Tyr1034 / Tyr1035)
Other Names	Phospho-Jak1 (Tyr1034 / Tyr1035)
Accession No.	Uniprot:P23458Gene ID:3716
Uniprot	P23458
GeneID	3716
Calculated MW	133
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only