

Phospho-MEK7 (Ser271 / Thr275) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91927Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-MEK7 (Ser271 / Thr275) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	VD(p-S)KAK(p-T)RS
Conjugates	Bio
Target Name	Phospho-MEK7 (Ser271 / Thr275)
Other Names	Phospho-MEK7 (Ser271 / Thr275)
Accession No.	Uniprot:O14733Gene ID:5609
Uniprot	O14733
GeneID	5609
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only