

Phospho-MLK3 (Thr277 / Ser281) Polyclonal Antibody HRP Conjugated



Catalog No: #C91928HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Phospho-MLK3 (Thr277 / Ser281) Polyclonal Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | K(p-T)TQM(p-S)AA |
| Conjugates | HRP |
| Target Name | Phospho-MLK3 (Thr277 / Ser281) |
| Other Names | Phospho-MLK3 (Thr277 / Ser281) |
| Accession No. | Uniprot:Q16584Gene ID:4296 |
| Uniprot | Q16584 |
| GeneID | 4296 |
| Calculated MW | 93 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only