

Phospho-MLK3 (Thr277 / Ser281) Polyclonal Antibody FITC Conjugated



Catalog No: #C91928FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-MLK3 (Thr277 / Ser281) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Ms Rat
Immunogen Type	K(p-T)TQM(p-S)AA
Conjugates	FITC
Target Name	Phospho-MLK3 (Thr277 / Ser281)
Other Names	Phospho-MLK3 (Thr277 / Ser281)
Accession No.	Uniprot:Q16584Gene ID:4296
Uniprot	Q16584
GeneID	4296
Calculated MW	93
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only