

Phospho-MYPT1 (Thr696) Polyclonal Antibody FITC Conjugated



Catalog No: #C91931FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-MYPT1 (Thr696) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	RS(p-T)QG
Conjugates	FITC
Target Name	Phospho-MYPT1 (Thr696)
Other Names	Phospho-MYPT1 (Thr696)
Accession No.	Uniprot:O14974Gene ID:4659
Uniprot	O14974
GeneID	4659
Calculated MW	113
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only