

phospho-SYN1 (Ser9) Polyclonal Antibody FITC Conjugated

Catalog No: #C91933FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-SYN1 (Ser9) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	RL(p-S)DS
Conjugates	FITC
Target Name	phospho-SYN1 (Ser9)
Other Names	phospho-SYN1 (Ser9)
Accession No.	Uniprot:P17600Gene ID:6853
Uniprot	P17600
GeneID	6853
Calculated MW	85
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only