

Phospho-PAK6 (Ser560) Polyclonal Antibody HRP Conjugated



Catalog No: #C91937HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-PAK6 (Ser560) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	RK(p-S)LV
Conjugates	HRP
Target Name	Phospho-PAK6 (Ser560)
Other Names	Phospho-PAK6 (Ser560)
Accession No.	Uniprot:Q9NQU5Gene ID:56924
Uniprot	Q9NQU5
GeneID	56924
Calculated MW	75
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000

Note: This product is for in vitro research use only