

Phospho-PDGFR α (Tyr1018) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91938Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-PDGFR α (Tyr1018) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	SG(p-Y)II
Conjugates	Bio
Target Name	Phospho-PDGFR α (Tyr1018)
Other Names	Phospho-PDGFR α (Tyr1018)
Accession No.	Uniprot:P16234Gene ID:5156
Uniprot	P16234
GeneID	5156
Calculated MW	117
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200?

Note: This product is for in vitro research use only