

Phospho-RPS6 (Ser240 / Ser244) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91951Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-RPS6 (Ser240 / Ser244) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	RA(p-S)TSK(p-S)ES
Conjugates	Bio
Target Name	Phospho-RPS6 (Ser240 / Ser244)
Other Names	Phospho-RPS6 (Ser240 / Ser244)
Accession No.	Uniprot:P62753Gene ID:6194
Uniprot	P62753
GeneID	6194
Calculated MW	29
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only