

# phospho-STAT1 (Ser727) Polyclonal Antibody APC Conjugated



Catalog No: #C91961APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | phospho-STAT1 (Ser727) Polyclonal Antibody APC Conjugated                       |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | IF  |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | PM(p-S)PE   |
| Conjugates         | APC   |
| Target Name        | phospho-STAT1 (Ser727)  |
| Other Names        | phospho-STAT1 (Ser727)  |
| Accession No.      | Uniprot:P42224Gene ID:6772  |
| Uniprot            | P42224  |
| GeneID             | 6772  |
| Calculated MW      | 87  |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only