

# phospho-STAT3 (Ser727) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91963Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-STAT3 (Ser727) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	PM(p-S)PR
Conjugates	Bio
Target Name	phospho-STAT3 (Ser727)
Other Names	phospho-STAT3 (Ser727)
Accession No.	Uniprot:P40763Gene ID:6774
Uniprot	P40763
GeneID	6774
Calculated MW	88
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only