

Phospho-VASP (Ser239) Polyclonal Antibody HRP Conjugated



Catalog No: #C91975HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-VASP (Ser239) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	KV(p-S)KQ
Conjugates	HRP
Target Name	Phospho-VASP (Ser239)
Other Names	Phospho-VASP (Ser239)
Accession No.	Uniprot:P50552Gene ID:7408
Uniprot	P50552
GeneID	7408
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only