

phospho-VEGF receptor 2 (Tyr1059) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91976Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-VEGF receptor 2 (Tyr1059) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	PD(p-Y)VR
Conjugates	Bio
Target Name	phospho-VEGF receptor 2 (Tyr1059)
Other Names	phospho-VEGF receptor 2 (Tyr1059)
Accession No.	Uniprot:P35968Gene ID:3791
Uniprot	P35968
GeneID	3791
Calculated MW	147
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only