

Phospho-mTOR (Ser2448) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91983Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospho-mTOR (Ser2448) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	TD(p-S)YS
Conjugates	Bio
Target Name	Phospho-mTOR (Ser2448)
Other Names	Phospho-mTOR (Ser2448)
Accession No.	Uniprot:P42345Gene ID:2475
Uniprot	P42345
GeneID	2475
Calculated MW	289
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only