

## FCGR1A/CD64 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91988Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FCGR1A/CD64 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/374
Conjugates	Bio
Target Name	FCGR1A/CD64
Other Names	FCGR1A/CD64
Accession No.	Uniprot:P12314Gene ID:2209
Uniprot	P12314
GeneID	2209
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only