

Caveolin-3 Polyclonal Antibody FITC Conjugated

Catalog No: #C92011FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caveolin-3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	95-151/151
Conjugates	FITC
Target Name	Caveolin-3
Other Names	Caveolin-3
Accession No.	Uniprot:P51637Gene ID:12391
Uniprot	P51637
GeneID	12391
Calculated MW	17
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only