

FPRL1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92021Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FPRL1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/351
Conjugates	Bio
Target Name	FPRL1
Other Names	FPRL1
Accession No.	Uniprot:P25090Gene ID:2358
Uniprot	P25090
GeneID	2358
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only