

## CD1d Polyclonal Antibody APC Conjugated

Catalog No: #C92033APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD1d Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Ms Rat
Immunogen Type	51-150/336
Conjugates	APC
Target Name	CD1d
Other Names	CD1d
Accession No.	Uniprot:P11609Gene ID:12479
Uniprot	P11609
GeneID	12479
Calculated MW	36
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only