

Neurogenin 2 Polyclonal Antibody FITC Conjugated

Catalog No: #C92034FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurogenin 2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/272
Conjugates	FITC
Target Name	Neurogenin 2
Other Names	Neurogenin 2
Accession No.	Uniprot:Q9H2A3Gene ID:63973
Uniprot	Q9H2A3
GeneID	63973
Calculated MW	30
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only