

Somatostatin Polyclonal Antibody HRP Conjugated

Catalog No: #C92038HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Somatostatin Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	103-116/116(disulfide 105,116)
Conjugates	HRP
Target Name	Somatostatin
Other Names	Somatostatin
Accession No.	Uniprot:P61278Gene ID:6750
Uniprot	P61278
GeneID	6750
Calculated MW	13
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only