

phospho-PKC delta (Thr505) Polyclonal Antibody FITC Conjugated

Catalog No: #C92041FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-PKC delta (Thr505) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	AS(p-T)FC
Conjugates	FITC
Target Name	phospho-PKC delta (Thr505)
Other Names	phospho-PKC delta (Thr505)
Accession No.	Uniprot:Q05655Gene ID:5580
Uniprot	Q05655
GeneID	5580
Calculated MW	77
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only