

ANO1/TMEM16A Polyclonal Antibody APC Conjugated

Catalog No: #C92056APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ANO1/TMEM16A Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/986
Conjugates	APC
Target Name	ANO1/TMEM16A
Other Names	ANO1/TMEM16A
Accession No.	Uniprot:Q5XXA6Gene ID:55107
Uniprot	Q5XXA6
GeneID	55107
Calculated MW	114
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only