

## GDF15 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92061Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | GDF15 Polyclonal Antibody Biotin Conjugated                                     |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | WB IHC IF   |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 221-308/308   |
| Conjugates         | Bio   |
| Target Name        | GDF15   |
| Other Names        | GDF15   |
| Accession No.      | Uniprot:Q99988Gene ID:9518  |
| Uniprot            | Q99988  |
| GeneID             | 9518  |
| Calculated MW      | 12.5/34   |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only