

Amphiregulin Polyclonal Antibody PE Conjugated

Catalog No: #C92065PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Amphiregulin Polyclonal Antibody PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IF |
| Species Reactivity | Hu Ms Rat |
| Conjugates | PE |
| Target Name | Amphiregulin |
| Other Names | Amphiregulin |
| Accession No. | Uniprot:P15514Gene ID:374 |
| Uniprot | P15514 |
| GeneID | 374 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF 1:100-500

Note: This product is for in vitro research use only