

ROR Gamma Polyclonal Antibody HRP Conjugated

Catalog No: #C92068HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ROR Gamma Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB?IHC
Species Reactivity	Hu Ms Rat
Immunogen Type	81-140/518
Conjugates	HRP
Target Name	ROR Gamma
Other Names	ROR Gamma
Accession No.	Uniprot:P51449Gene ID:6097
Uniprot	P51449
GeneID	6097
Calculated MW	57
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only