

Cytochrome b245 Light Chain Polyclonal Antibody PE Conjugated



Catalog No: #C92072PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytochrome b245 Light Chain Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/195
Conjugates	PE
Target Name	Cytochrome b245 Light Chain
Other Names	Cytochrome b245 Light Chain
Accession No.	Uniprot:P13498Gene ID:1535
Uniprot	P13498
GeneID	1535
Calculated MW	21
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only