

ADCY7 Polyclonal Antibody APC Conjugated

Catalog No: #C92088APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ADCY7 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	901-1080/1080
Conjugates	APC
Target Name	ADCY7
Other Names	ADCY7
Accession No.	Uniprot:P51828Gene ID:113
Uniprot	P51828
GeneID	113
Calculated MW	120
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only