

phospho-AMPK alpha-2 (Thr172) Polyclonal Antibody Biotin Conjugated



Catalog No: #C92102Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-AMPK alpha-2 (Thr172) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	LR(p-T)SC
Conjugates	Bio
Target Name	phospho-AMPK alpha-2 (Thr172)
Other Names	phospho-AMPK alpha-2 (Thr172)
Accession No.	Uniprot:P54646Gene ID:5563
Uniprot	P54646
GeneID	5563
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only