

DRP1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92121Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DRP1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/763
Conjugates	Bio
Target Name	DRP1
Other Names	DRP1
Accession No.	Uniprot:O00429Gene ID:10059
Uniprot	O00429
GeneID	10059
Calculated MW	81
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only