

MST1R Polyclonal Antibody FITC Conjugated

Catalog No: #C92130FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MST1R Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/1400
Conjugates	FITC
Target Name	MST1R
Other Names	MST1R
Accession No.	Uniprot:Q04912Gene ID:4486
Uniprot	Q04912
GeneID	4486
Calculated MW	30/119/150
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only