

MST1R Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C92130PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | MST1R Polyclonal Antibody PE-Cy7 Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 151-250/1400 |
| Conjugates | PE-Cy7 |
| Target Name | MST1R |
| Other Names | MST1R |
| Accession No. | Uniprot:Q04912Gene ID:4486 |
| Uniprot | Q04912 |
| GeneID | 4486 |
| Calculated MW | 30/119/150 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only