

VAV1 Polyclonal Antibody FITC Conjugated

Catalog No: #C92138FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VAV1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	741-845/845
Conjugates	FITC
Target Name	VAV1
Other Names	VAV1
Accession No.	Uniprot:P15498Gene ID:7409
Uniprot	P15498
GeneID	7409
Calculated MW	98
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only