

## CD167b/DDR2 Polyclonal Antibody PE Conjugated

Catalog No: #C92148PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD167b/DDR2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	245-350/855
Conjugates	PE
Target Name	CD167b/DDR2
Other Names	CD167b/DDR2
Accession No.	Uniprot:Q16832Gene ID:4921
Uniprot	Q16832
GeneID	4921
Calculated MW	92
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only