

Neurologin 1 Polyclonal Antibody APC Conjugated

Catalog No: #C92151APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurologin 1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	APC
Target Name	Neurologin 1
Other Names	Neurologin 1
Accession No.	Uniprot:Q8N2Q7Gene ID:22871
Uniprot	Q8N2Q7
GeneID	22871
Calculated MW	87
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only